

BD Molecular Swab Collection Kit: Vaginal swab specimen clinician collection and transfer procedure


Clinician collection procedure

1. Collect swab prior to pelvic, speculum, or bimanual exam.
2. Gently slide the swab no more than 2 inches (5 cm) into the vagina. Do not use lubricants or other products containing substances such as carbomers.
3. Rotate the swab for 10 to 15 seconds.
4. Withdraw the swab without touching the skin outside the vagina.


Swab-to-tube transfer procedure

Specimens collected using the BD Molecular Collection Swab must be transferred to the BD Molecular Swab Sample Buffer Tube immediately after collection.


To transfer the sample




1. Unscrew the cap of the BD Molecular Swab Sample Buffer Tube, taking care not to contaminate the contents or the outside of the tube. Immediately after collection, insert the BD Molecular Collection Swab into the tube so that the score mark indicated by the black line is at the lip of the tube.



2. Carefully break the shaft at the score mark and allow the swab to drop into the tube.



3. Tightly re-cap the tube.



4. Label tube with patient information, date, and time collected. Be careful not to obscure the barcodes on the tube.

Assay	Condition	Duration
BD CTGC2 for BD MAX™ System BD CTGCTV2 for BD MAX™ System	2 – 30 °C	Up to 21 days
BD MAX™ Vaginal Panel	2 – 8 °C	Up to 14 days
	2 – 30 °C	Up to 8 days

Approved for use with:

- BD CTGC2 for BD MAX™ System
- BD CTGCTV2 for BD MAX™ System
- BD MAX™ Vaginal Panel

BD - Europe, Terre-Bonne Park - A4, Route de Crassier 17, 1262 Eysins, Switzerland

The BD MAX™ System, BD Vaginal Panel and BD Molecular Swab Collection Kit are in vitro diagnostic medical devices bearing a CE-mark. The BD CTGCTV2 & BD CTGC2 for BD MAX™ System are in vitro diagnostic medical devices bearing a CE mark and is CE certified by BSI Group The Netherlands B.V. (Notified Body Number = 2797)

bd.com



BD, the BD Logo and MAX are trademarks of Becton, Dickinson and Company or its affiliates. © 2021 BD. All rights reserved. BD-39665

